

SANO₃

OZONE THERAPY KIT FOR
MAJOR AUTOHEMO INFUSION
Patented

MADE IN A PHARMACEUTICAL PLANT



CE Marked

OZONE THERAPY KIT FOR MAJOR AUTOHEMO INFUSION

CLOSED CIRCUIT

SANO3 is a closed circuit system that allows all Major Autohemo Infusion procedures to be performed under extremely safe conditions, avoiding risks of environmental contamination. This safety also extends to the health care provider and the patient.

PHTHALATE-FREE MATERIAL

All materials constituting the SANO3 system are phthalate-free and certified for 100% compatibility with the O_2O_3 mixture.

NO NEEDLE

The SANO3 system stands apart as the needle is only required for the vein connection. The needle-free system makes it easy to insert the O_2O_3 mixture into the bag. Indeed, the use of a syringe without a needle eliminates the risk of accidental puncture. This approach ensures a safer experience for medical staff and patients, improving the efficiency of medical practice.

AUTOMATIC CONNECTIONS

The SANO3 system is equipped with advanced proprietary needle-free on-off connections for tight sealing both before and after use.

ABSENCE OF BLOOD CLOTS AND GAS EMBOLISM


SANO3 allows Major Autohemo Infusion procedures to be performed without the risk of blood clots and gas embolism.

SANO₃

THE KIT

Each Major Autohemo Infusion Kit consists of two medical devices individually marked CE.

1. Sterile 400 mL PVC bag, DEHP-free, containing 35 mL of ACD-A anticoagulant solution.
2. Sterile transfusion set with 175-micron filter, equipped with Y-injection port and Luer Lock connection.



For proper use, follow the directions in the instructions for use.



SALES DEPARTMENT
I - 20135 Milano (MI)
Viale Lazio, 26
Tel +39 02 5513513
Web: haemopharm.it



HEAD OFFICE and PRODUCTION FACILITY
I - 23030 Tovo di Sant'Agata (SO)
Via Provinciale, 8
Tel +39 0342 771017
E-mail: info@haemopharm.it

All information in this presentation and about the product described may be changed at any time by Haemopharm Healthcare without prior notice.
Information regarding the specific product in use is only that included in the product packaging.

SANO3 is a patented system composed of two individually CE marked medical devices produced by Paolo Gobbi Frattini S.r.l.
SANO3 is a system distributed by Haemopharm Healthcare S.r.l.

©Copyright 2024 Haemopharm Healthcare S.r.l. All rights reserved. REV. 03/ 05_2024